


magnetic charger kit




off



low power



full power



i-limb® magnetic battery charger kit with LED battery display

- Simple to use on/off switch
- Low battery LED warning light
- Fail-safe feature - battery reserve prevents hand locking in closed position

ÖSSUR

TOUCH
SOLUTIONS

Össur Touch Solutions are designed to ensure the best possible outcomes for people with upper-limb deficiencies.



WWW.OSSUR.COM

Össur Europe B.V.
De Schakel 70
5651 GH Eindhoven
P.O. Box 120
5690 AC Son en Breugel
The Netherlands

TEL +800 3539 3668
+31 (0) 499 462 840
FAX +800 3539 3299
+31 (0) 499 462 841
info-europe@ossur.com

 **ÖSSUR**
LIFE WITHOUT LIMITATIONS