

# flexion wrist

Enjoy the freedom and flexibility of natural wrist movement



## Features:


- Smooth & natural motion which provides control in completion of daily tasks
- Enables shock and torque absorption
- Lock and unlock function via simple magnetic lock
- Particularly beneficial for bilateral patients
- Available in both black and neutral colors to match your **i-limb** hand



Lock **flexion wrist** to position **i-limb** at  $+40^\circ$ ,  $+20^\circ$ ,  $0^\circ$ ,  $-20^\circ$ ,  $-40^\circ$



Unlock **flexion wrist** for spring-loaded flexibility range from  $+40^\circ$  to  $-40^\circ$



Locking **flexion wrist** in various positions allows for better positioning during functional tasks such as:


- Writing
- Brushing hair
- Typing
- Using utensils for cooking and eating
- Fastening clothes
- Manipulating small objects

Unlocking **flexion wrist** allows for better positioning during functional tasks such as:

- Sweeping
- Cleaning tabletops
- Pushing shopping cart
- Opening doors
- Stabilizing objects during cooking tasks

# flexion wrist

All lengths stated in mm



Medium i-limb with flexion wrist

Small i-limb with flexion wrist

## Technical Information

flexion wrist weight	161g
Range of movement (unlocked)	+40° to -40°
Locked positions	+40°, +20°, 0°, -20°, -40°

Product Number	Standard Hand
PL270001	i-limb ultra revolution flexion wrist Medium RH Neutral
PL271001	i-limb ultra revolution flexion wrist Medium LH Neutral
PL272001	i-limb ultra revolution flexion wrist Small RH Neutral
PL273001	i-limb ultra revolution flexion wrist Small LH Neutral
PL274001	i-limb ultra revolution flexion wrist Medium RH Black
PL275001	i-limb ultra revolution flexion wrist Medium LH Black
PL276001	i-limb ultra revolution flexion wrist Small RH Black
PL277001	i-limb ultra revolution flexion wrist Small LH Black

Product Number	Standard Hand
PL206001	i-limb ultra flexion wrist Medium RH Neutral
PL207001	i-limb ultra flexion wrist Medium LH Neutral
PL208001	i-limb ultra flexion wrist Small RH Neutral
PL209001	i-limb ultra flexion wrist Small LH Neutral
PL210001	i-limb ultra flexion wrist Medium RH Black
PL211001	i-limb ultra flexion wrist Medium LH Black
PL212001	i-limb ultra flexion wrist Small RH Black
PL213001	i-limb ultra flexion wrist Small LH Black

North American Customers (US, Canada & Mexico)  
Tel: +1 855 MY iLIMB (694 5462)

Europe & International Customers  
Tel: +44 1506 438 556

www.touchbionics.com  
Email: info@touchbionics.com