


precision. power. intelligent motion.




## Introducing i-limb quantum:

- **Smarter:**  
i-mo™ technology uses simple gestures to change grips
- **Faster:**  
Boost digit speed by up to 30%
- **Stronger:**  
Up to 30% more power when needed
- **Smaller:**  
Anatomical styling now available in 3 sizes
  
- **Speed boost:**  
Adjustable open/close speed for digits



All lengths stated in mm



Technical Information	
Voltage	7.4 V (nominal)
Max. Current	5 A
Battery Capacity	Rechargeable lithium polymer; 7.4 V (nominal); 2000 mAh capacity; 1300 mAh capacity
Max. hand load limit (static limit)	40kg/88lb (Extra Small) 90kg/198lb (Small/Medium)
Finger carry load (static limit)	20kg/44lbs (Extra Small) 32kg/71lbs (Small/Medium)
Time from open position to full power grip	0.8 seconds
Device Weight	Weight with QWD 454g/1lb (Extra Small); 504g/1.11lbs (Small); 516g/1.14lbs (Medium) Weight with WD 418g/0.92lbs (Extra Small); 467g/1.03 lbs (Small); 479g /1.06lbs (Medium) Weight with flexion wrist 561g/1.24lbs (Extra Small); 611g/1.35lbs (Small); 623g/1.37lbs (Medium) Weight with friction wrist 449g/0.99lbs (Extra Small); 499g/1.10 lbs (Small); 511g/ 1.13lbs (Medium)

North American Customers (US, Canada & Mexico)  
Tel: +1 855 MY iLIMB (694 5462)

Europe & International Customers  
Tel: +44 1506 438 556

[www.touchbionics.com](http://www.touchbionics.com)  
Email: [info@touchbionics.com](mailto:info@touchbionics.com)