

Coverings are an important part of the appearance and durability of upper limb prostheses. Touch Bionics i-limb products are the first prosthetic hands and fingers to imitate the individual digit articulation of the human hand, and we provide flexible coverings which enable that articulation. Our covering solutions include:

## i-limb<sup>tm</sup> skin active

For users who prefer the robotic nature of the uncovered i-limb prosthesis, **i-limb skin active** is a thin layer of semi-transparent or black material that has been designed to conform to every contour of the hand.



## i-limb<sup>tm</sup> skin active TS

In addition to all the benefits of the **i-limb skin active**, the **i-limb skin active TS** allows users to use touchscreen devices such as smartphones or tablets with compatible i-limb prostheses.



## i-limb<sup>tm</sup> skin natural

Providing the user a lifelike silicone covering, **i-limb skin natural** is available with 18 different skin shades.

Benefits include:

- Don or doff in under a minute
- Improved design optimizes hand power and precision pinch mode



## i-limb<sup>tm</sup> skin match

For the most realistic covering option possible, this is a fully customized silicone covering that exactly matches the user's skin tone and features.

