

livingskin

Passive functional prostheses made with high definition silicone.

- Custom sculpting and painting ensure an individualized fit including:
 - Optimum match to patient skin tone
 - Personalized detail with hair, freckles, moles or scarring
- Functional capabilities can assist in Activities of Daily Living such as pushing, pulling, stabilizing, supporting, light grasping, and typing
- Skin colour and other features are beneath the outer clear silicone layer to prevent fading
- SeasonGuard® can help adjust the color to ensure match throughout the year
- Adjustable armatures may be fitted to allow repositioning of each finger for specific tasks
- Silicone provides a versatile, durable and virtually stain-proof prosthetic option
- Complementary solution to a myoelectric prosthesis to alternate with passive prosthesis as desired



various options available



ÖSSUR

TOUCH SOLUTIONS

Össur Touch Solutions are designed to ensure the best possible outcomes for people with upper-limb deficiencies.

ACCEPTANCE GUARANTEE



Touch Bionics is committed to deliver the highest quality livingskin products. Our Acceptance Guarantee ensures that the product is delivered to meet your patient's needs.

ORDERING A LIVINGSKIN PROSTHESIS

We offer the following program for clinicians who would like to fit a livingskin prosthesis.

MY LIVINGSKIN

Our revolutionary iPad-based color capture system enables you to accurately match skin tones and select features such as hair, nails, veins, or freckles. With an intuitive design, the easy on screen step-by-step instruction guides you through the entire process resulting in a high quality prosthesis or customized i-limb skin match.



FOLLOW US ON



USA (855) 694-5462
CANADA (416) 465-0266
WWW.TOUCHBIONICS.COM

