

virtu-limb is a wireless simulation and training product for myoelectric upper limb prostheses.

Advanced connectivity:

- Connects wirelessly to your computer via Bluetooth software
- Mobile app (available through the Apple® app store) allows for remote EMG testing in any location

Features:

- **virtu-limb** can be used with **i-limb** products and other myoelectric devices
- Option of using an on-screen virtual representation of an **i-limb ultra revolution**
- Compatible with compact electrodes, remote electrodes or FSR touchpad
- Compatible with quick wrist disconnect and **flexion wrist** options

virtu-limb can be used for a variety of purposes:

- Measuring patients' myoelectric signals to determine the best muscle sites for their prosthesis
- Training patients how to control their muscle signals
- Training patients how to use features
- Demonstrating **i-limb** products to patients

